

FIAlab Instruments Inc.



Fluidics Intelligently Automated



HELLO!

I am Scott Schroeder
National Sales Manager

You can find me at scott@flowinjection.com
Home Office (262-670-1389) Cell (206-619-6705)



THEY LOVE FLOW INJECTION



Fluidics Intelligently Automated

FIAlab Instruments has perhaps the strongest roots in flow injection out there. The inventor of these techniques (J. Ruzicka) founded FIAlab in 1987. For over 30 years now, we have remained on the cutting-edge of these techniques and have significantly improved these technologies.

FIAlab Improvements over Continuous Flow Analyzers

- Size-Bench Space Critical
- Backup for Lab (redundancy)
- Full Array Spectrometer
- Pump Control (FIA)
- Method Improvements
(Hardware/Software/Chemistry)
- Significant Reduction in cost of operation

EXTREMELY SMALL SIZE 25 cm X 15.2 cm



60 cm X 25 cm X 50 cm
Reagents 5 cm (actual space 10cm)

SIMULTANEOUSLY RUN MULTIPLE ANALYTES =8



Ammonia/TKN all on 1 Channel no switchover either FIA100 or EPA 350.1/351.2
NO₃+NO₂ 353.2

UNIQUE FLEXIBILITY WITH LITTLE ADDED COST



■ FLEXIBILITY

■ Add Autosampler and Software license SPLIT

■ TWO WORKSTATIONS

■ THREE WORKSTATIONS

HIGH END SPECTROMETER IN FIA1000

Flame Spectrometer


The Flame Spectrometer is built using industry-leading manufacturing techniques that help deliver high thermal stability and low unit to unit variation- without compromising the flexibility and configurability that are the hallmark of Ocean Insight's miniature spectrometers. New features such as interchangeable slits, indicator LEDs, and simple device connectors deliver more freedom and less frustration.

Whether you are a research lab looking to make a breakthrough, or an engineer working to integrate a spectrometer into an OEM system, the Flame will provide you with the performance and features you need to make your UV-Vis spectroscopy application successful.




ALL WAVELENGTHS SIMULTANEOUSLY

Ammonia



NO3+NO2



SOFTWARE CONTROLLED ISMATEC PUMP



FL-1000 Pump

The FL-1000 pump features a planetary drive for less pulsation and longer tubing life. The eight SS rotors provides smooth rotation and accurate liquid flow. This pump requires fewer tubing changes and can handle flow rates of 0.002 to 44 mL/min per channel. Flow rates depend on drive rpm and tubing size. See page XX for a complete list of FIAlab's pump tubing. The FL-1000 pump is simple to use and includes automatic occlusion cartridges, which provide reproducible results.

Pump Extensions

The FL-1000 pump can be extended up to 8 channels to accommodate multiple chemistries.



Fluorometric Ammonia & TKN


Alkalinity

Method number	Lower	Upper	MDL	Units	Sample / Hour	Matrix	Compliant With	Notes
ALK-W-1-1	1 10	50 500	0.3 3	g CaCO ₃ / L	50	Waters	EPA 310.2	Methyl Orange method.

Ammonia

Method number	Lower	Upper	MDL	Units	Sample / Hour	Matrix	Compliant With	Notes
NH3-S-1-1	0.02 0.075	0.5 50	0.006 0.025	mg N / L as NH ₃	120	Soil extracts	N/A	Salicylate method for soil extracts.
NH3-W-1-2	0.5	20	0.1	mg N / L as NH ₃	60	Waters	EPA 350.1	Salicylate method with gas diffusion.
NH3-W-1-4	0.5	20	0.1	mg N / L as NH ₃	60	Waters	SM 4500-NH ₃ H.	Salicylate or phenate method with gas diffusion.
NH3-W-2-1	0.05	0.5	0.02	mg N / L as NH ₃	60	Waters	EPA 350.1	Salicylate or phenate method with gas diffusion, utilizing low-noise detector.
NH3-W-2-3	0.05	0.5	0.02	mg N / L as NH ₃	60	Waters	SM 4500-NH ₃ H.	Salicylate or phenate method with gas diffusion, utilizing low-noise detector.
NH3-W-2-5	0.006	1	0.002	mg N / L as NH ₃	60	Waters	SM 4500-NH ₃ H.	Salicylate or phenate method, utilizing low-noise detector.
NH3-W-3-2	0.05	10	0.012	mg N / L as NH ₃	60	Waters	EPA Rec.	OPA method with gas diffusion, utilizing fluorometric detector.

•APPROVED
JULY 19th 2021



Gas Diffusion works well Ammonia



FIALAB VERSION

LESS REAGENT CONSUMPTION & WASTE GENERATION

Chemistry	Lachat method	Standard Method	Fialab method	Low	High	units	Lachat	FIALab
Nitrate/Nitrite	10-107-04-1-C	4500-N03-I	NO3-W-20-2	0.2	2	mg/L		
Nitrate/Nitrite		4500-N03-I		0.2	2	mg/L		
Chemical	Chem Formula		Reagent	Reagent Name		prepare		
Ammonium chloride	NH ₄ Cl		1	1.6 M Ammonium Chloride Buffer (1 L)		3 months	3.07	1.925
Sodium Hydroxide	NaOH		2	Sulfanilamide Solution (1 L)		Monthly	0.94	0.94
Ethylenediaminetetraacetic acid disodium salt dihydrate	C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ · 2H ₂ O							
Sulfanilamide	H ₂ NC ₆ H ₄ SO ₂ NH ₂		3	Carrier Water or 2 ml H ₂ SO ₄ /L/4ml FIA			1.28	0.77
Phosphoric acid, 85%	H ₃ PO ₄			Total (ml)/minute			5.29	3.635
N-1-Naphthylethylene diamine dihydrochloride	C ₁₀ H ₇ NHCH ₂ CH ₂ NH ₂ · 2HCl							
Brij® L23, 30% solution	CH ₃ (CH ₂) ₁₀ CH ₂ (OCH ₂ CH ₂) _n OH			Waste comparison Reduction of waste approximately 30%				
Sulfuric acid (36N)	H ₂ SO ₄							

FLOW CELLS

STS Series Spectrometers
Ocean Optics

Flow Cells

Professionally machined. Carefully measured. Proudly made in the USA. These flow cells come in a wide range of materials and lengths to satisfy your detection limits and withstand your reagents. Normal and micro-volumes are available.

Products Include:

- SMA-Z Flow Cells 100um to 100mm optical length
- Long Path 50 cm optical path
- Dialysis and gas diffusion cells




Assorted flow cells

Materials Available:

- Plexiglas
- Ultem
- Teflon
- Stainless Steel
- Peek
- COP

LAB ON VALVE DESIGN



DEMO TIME



Fluidics Intelligently Automated



Conclusions



THANK YOU!